


Assign segmentation to article

To assign a segmentation to an article, follow these steps:

1. Click , the following window appears:



The screenshot shows a 'New Record' dialog box. The title bar is blue and contains the text 'New Record' and a close button (X). The main area is white and contains a text input field with the label 'Segmentation Id *'. At the bottom of the dialog, there are three buttons: 'CANCEL', 'SAVE', and 'SAVE AND NEW'.

You can also assign a segmentation by selecting “Create Segmentation” from the Actions menu next to article.

2. Click the “Segmentation Id” field then select the segmentation you want,
3. Click **Save** to add the segmentation and close the window, or **Save and New** to add it and open a new window. The segmentation is added to the list.

You can also assign a segmentation to multiple articles as mentioned in “[Set Segmentation to Articles](#)”.

Note:

- This button is not available if the selected catalog is a child catalog with the “Segmentation Mirroring” option activated.